

Scaling-up strategies for controllable biosynthetic ZnO NPs using cell free-extract of endophytic *Streptomyces albus*: characterization, statistical optimization, and biomedical activities evaluation

Shahira H. EL-Moslami, Mohamed S. Elnouby, Ahmed H. Rezk & Esmail M. El-Fakharany

Abstract

A cell-free extract of endophytic *Streptomyces albus* E56 was used to biosynthesize ZnO nanoparticles under fed-batch conditions. Process optimization using Plackett–Burman and response surface methods increased lab-scale yield from 4.36 to 7.59 g/L (~1.6-fold). Scale-up using a 7 L pre-culture and a 70 L stainless steel bioreactor with exponential pulse feeding achieved a 31.25-fold increase (~155 g/L). The resulting ZnO/MnO₂ composite showed strong antimicrobial activity.



Winpact Model: FS-07 & FS-V-D

Introduction

Microbial biosynthesis of ZnO nanoparticles offers eco-friendly and controllable alternatives to chemical methods. However, there's a gap in scalable processes. This study fills it by using fed-batch fermentation for bioprocess scale-up from lab to semi-industrial level with characterized antimicrobial efficacy..

Materials and Methods

Cultivation was firstly optimized in lab fed-batch mode. Scale-up involved a stirred 7 L pre-culture followed by fed-batch fermentation in a 70 L stainless steel bioreactor, applying exponential pulse feeding with a 1:200 N/C ratio over 120–312 h. Dissolved oxygen was maintained above 20% via airflow and agitation.

US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
Taoyuan City 330041, Taiwan
T/+886-3-3762878
F/+886-3-3761310

Shanghai Office

Room 2519, No. 306 Taopu Road,
Putuo District, Shanghai, China
Postal code: 200333
Cell:+86-13621808204
T/ +86-21-50795277
F/ +86-21-50795277

India Office

D.No.7-143, 2nd Floor,
St.No.2,Nagendra Nagar,
Habsiguda, Hyderabad-500007,
India
T/ +91-40-27001515
T/ +91-40-27001586



www.majorsci.com
info@majorsci.com

Results

- Lab-scale yield increased $\sim 1.6\times$ ($4.36 \rightarrow 7.59$ g/L ZnO NPs).
- Pilot scale fed-batch in 70 L reactor increased yield $31.25\times$ (to ~ 155 g/L).
- SEM/TEM showed 6–18 nm spherical particles.
- Antimicrobial testing revealed substantial inhibition zones against drug-resistant bacteria and pathogens (e.g., *Salmonella paratyphi* inhibition zone ~ 53 mm).

References

Scaling-up strategies for controllable biosynthetic ZnO NPs using cell-free extract of endophytic *Streptomyces albus*: characterization, statistical optimization, and biomedical activities evaluation

<https://www.nature.com/articles/s41598-023-29757-9>

US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
Taoyuan City 330041, Taiwan
T/+886-3-3762878
F/+886-3-3761310

Shanghai Office

Room 2519, No. 306 Taopu Road,
Putuo District, Shanghai, China
Postal code: 200333
Cell: +86-13621808204
T/ +86-21-50795277
F/ +86-21-50795277

India Office

D.No.7-143, 2nd Floor,
St.No.2,Nagendra Nagar,
Habsiguda, Hyderabad-500007,
India
T/ +91-40-27001515
T/ +91-40-27001586