

3D Mini Tray Shaker

MS-3D-TS

 major science

www.majorsci.com
info@majorsci.com

產品影片





NEW



Specification

Name	3D Mini Tray Shaker
Cat. No.	MS-3D-TS
Angle	20 °
Movement	Three-Dimensional Movement
Speed range	0~27rpm
Speed display	scale (min/max)
Platform size(WxL)	Approx. 9.21"x6.61"(168x234mm)
Overall size(WxLxH)	Approx. 6.93"x8.66"x6.5"(176×220×165mm)
Max Loading	2 kg
Weight	1.7 kg
Power	20W
Input voltage	AC100~240V 50/60Hz



MS-3D-TS-001

Optional Accessory

Cat. No.	Product Description
MS-3D-TS-001	3D Spring Tray

* All images are for reference only, actual products might differ from the pictures above.
* Technical specifications subject to change without notice.

Innovative Life Sciences Tools

Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
Taoyuan City 330041, Taiwan
T/+886-3-3762878
F/+886-3-3761310

US Office


19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

Shanghai Office

Room 612, International business exhibition center,
9300 Hunan Road, Pudong, Shanghai, China
National toll-free No.:400-823-9177
T/ +86-21-50795277
F/ +86-21-50795277

India Office

D.No.7-143, 2nd Floor,
St.No.2,Nagendra Nagar,
Habsiguda, Hyderabad-500007.
India
T/ +91-40-27001515
T/ +91-40-27001586



*Could be compatible with MO-MINI.